


COMBI TOUCH *Spare Parts*



Spare Parts Price List

Nr.	Description	Traceability	Warranty	Price
00400001	IRRIGATION TUBE - 0,8 mt each unit		CONSUMABLE	€ 9,80 /mt
00400007	AIR OUTPUT SIATPUR TUBE - 0,7 mt each unit		CONSUMABLE	€ 19,90 /mt
00400014	AIR INPUT TUBE - 1 mt each unit		CONSUMABLE	€ 4,80 /mt
00410014	O-RING 2,90X1,78 VITON		CONSUMABLE	€ 2,10
02950005-001	O-RING 2,90X1,78 VITON - KIT 20Pcs		CONSUMABLE	€ 20,00
01640008	MICROSWITCH OMRON D2FL		6 MONTHS	€ 8,64
01640010	ON-OFF GREEN SWITCH		6 MONTHS	€ 19,20
01980001	RUBBER FOOT		6 MONTHS	€ 2,00
02950026-001	RUBBER FOOT - KIT 200Pcs		6 MONTHS	€ 60,00
02790115	ACTUATOR REED		6 MONTHS	€ 12,00
02820065-001	COMPLETE BASE BOARD WITH CABLE WIRING	LOT	6 MONTHS	€ 320,00
02830055-001	LED SCALER CORD	LOT	6 MONTHS	€ 210,00
03020050-001	WATER TAP - RIGHT HAND SIDE		6 MONTHS	€ 53,00
03020193-001	3-WAY VALVE		6 MONTHS	€ 145,00
03020196-001 ³	COMPLETE HEATER		6 MONTHS	€ 410,00
03020217-001	WATER TAP - LEFT HAND SIDE		6 MONTHS	€ 82,00
03020218-001 ²	MOTOR	LOT	6 MONTHS	€ 85,00
03020223-001	SCALER HANDPIECE SUPPORT		6 MONTHS	€ 54,30
03020224-001	JET POLISHER HANDPIECE SUPPORT		6 MONTHS	€ 61,00
03020226-001	PERISTALTIC PUMP BLACK		6 MONTHS	€ 79,00
03020242-001	TUBES 2X4 Y-SHAPED COUPLING		6 MONTHS	€ 64,90
03020246-001	FAN		6 MONTHS	€ 52,00
03020267-001	POWDER CONTAINER AIR INLET KIT		6 MONTHS	€ 63,20
03020268-001	POWDER CONTAINER AIR OUTLET KIT		6 MONTHS	€ 73,50
03020289-001 ³	COMPLETE AIR GROUP	LOT	6 MONTHS	€ 900,00
03020290-001 ³	ASSEMBL. AIR INPUT SOLENOID VALVE		6 MONTHS	€ 78,40
03070069-001	COMPLETE COVER WITH TOUCH SCREEN	LOT	6 MONTHS	€ 620,00
03100018-001	CONTROL MODULE	LOT	6 MONTHS	€ 420,00
03120222-001	JET POLISHER HANDPIECE WITH CORD		6 MONTHS	€ 290,00
03150103 ¹	WHITE PUMP SUPPORT		6 MONTHS	€ 10,00
03210001-001	PERISTALTIC PUMP WHITE		6 MONTHS	€ 79,00
04430002-001	COMPLETE POWER SUPPLY	LOT	6 MONTHS	€ 270,00



Notes

- (1) In case you have a unit with a BLACK peristaltic pump and you change it to the WHITE one, you will need to change also the pump support and order the item 03150103.
- (2) The item 03020218-001 includes motor and plastic support for the WHITE pump. In case you have a unit with a BLACK peristaltic pump you will not need to use the plastic support supplied.
- (3) Refer to Service Notes nr 2020/001 and 2020/002 for the upgrades on these items